



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Rabbit Anti-Murine PAI-1 Anti-Serum

Product #: PI 35

Lot #: ASMPAI-896M

Molecular Weight: 160 KDa

Physical Appearance: colorless frozen solution

Solubility: <10 mg/mL

Concentration: 10 mg/mL

Applications:

ELISA: 1:10,000 to 1:100,000 dilution

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr

Buffer: 0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 6.6