



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse PAI-1 Fluorescent Labeled IgG
Product Number: PI 38
Aliquot: 1.0 mg
Lot Number pi38.1105A
Storage: -80°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is specific for mouse PAI-1. This is a fluorescent labeled IgG purified fraction.
MOLECULAR WEIGHT:	43 kDa non-glycosylated PAI-1 50 kDa glycosylated PAI-1
CONCENTRATION:	2.8 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS. PROTECT FROM LIGHT.
IMMUNOGEN:	Recombinant Mouse PAI-1
ISOTYPE:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
PURITY:	Not Determined.
APPLICATIONS:	Western blot: 1-5 µg/mL Immunohistochemistry: 1-10 µg/mL
SPECIES SPECIFICITY:	Mouse.