



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Rat PAI-1 Antiserum

Product Number: PI 39

Aliquot: 1.0 mL

Lot Number: pi39.103A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for rat plasminogen activator inhibitor-1 (PAI-1) and is provided as neat antiserum.
MOLECULAR WEIGHT:	50 kDa glycosylated 43 kDa non-glycosylated
CONCENTRATION:	Not determined
STORAGE BUFFER:	Not applicable
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Rat PAI-1
PURIFICATION:	Not purified
APPLICATIONS:	Suitable for use in Western blotting, ELISA and IHC applications.
SPECIES SPECIFICITY:	Rat