



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Rat PAI-1 IgG Fraction

Product Number: PI 40

Aliquot: 1.0 mg

Lot Number: pi40.1299A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for rat plasminogen activator inhibitor-1 (PAI-1) and is provided as purified IgG.
MOLECULAR WEIGHT:	50 kDa glycosylated 43 kDa non-glycosylated
CONCENTRATION:	10.0 mg/mL
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.15 M NaCl; pH 7.0
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Rat PAI-1
PURIFICATION:	Purified by Protein A Sepharose
APPLICATIONS:	Suitable for use in Western blotting, ELISA and IHC applications.
SPECIES SPECIFICITY:	Rat