



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Rat PAI-1 IgG – Biotin Labeled

Product Number: PI 41

Aliquot: 1.0 mg

Lot Number: pi41.103A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for rat plasminogen activator inhibitor-1 (PAI-1). This IgG prep has been biotin labeled for greater utility.
MOLECULAR WEIGHT:	50 kDa glycosylated 43 kDa non-glycosylated
CONCENTRATION:	7.84 mg/mL as determined spectrophotometrically.
STORAGE BUFFER:	0.1M HEPES; 0.1M NaCl; 1mM EDTA; pH 7.4
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Rat PAI-1
PURIFICATION:	Protein A Sepharose purification
CLASS:	IgG
APPLICATIONS:	Suitable for use in Western blotting, ELISA and IHC applications.
SPECIES SPECIFICITY:	Rat