



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Rat PAI-1 Affinity Purified IgG

Product Number: PI 43

Aliquot: 100 µg

Lot Number: pi43.103A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for rat plasminogen activator inhibitor-1 (PAI-1) and has been affinity purified for improved quality.
MOLECULAR WEIGHT:	50 kDa glycosylated 43 kDa non-glycosylated
CONCENTRATION:	10.98 mg/mL as determined spectrophotometrically.
STORAGE BUFFER:	0.1M HEPES; 0.1M NaCl; 1mM EDTA; pH 7.4
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Rat PAI-1
CLASS:	IgG
PURIFICATION:	Affinity purified by immobilized PAI-1.
APPLICATIONS:	Suitable for use in Western blotting, ELISA and IHC applications.
SPECIES SPECIFICITY:	Rat