



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

PAI-1 Depleted Rat Plasma (in Na Citrate)

Product Number: PI 52

Aliquot: 10 mL

Lot Number pi52.306A

Storage: -80°C

DESCRIPTION: This rat plasma preparation is depleted of PAI-1 by using immobilized anti-rat PAI-1 antibody.

**ACTIVE PAI-1
CONCENTRATION:** 0.07 ng/mL

PRESENTATION: Rat plasma in Na Citrate void of PAI-1.

STORAGE: -80°C