



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Sheep Anti-Human PAI-1 Affinity Purified
Product Number: PI59
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Anti-PAI-1 (Human), Sheep Polyclonal, Affinity Purified
MOLECULAR WEIGHT:	160,000
CONCENTRATION:	10 mg/mL
SIZE:	0.1 mg
VOLUME:	0.01 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Sheep serum
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
PURIFICATION:	Affinity purified by immobilized human PAI-1

For research purposes only.