



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Mouse Anti-Human PAI-1
Product Number: PI62
Aliquot: 1.0 mg
Lot Number: pi62.1209
Storage: -70°C

DESCRIPTION:	This is an inhibitory monoclonal antibody to Human PAI-1. It also cross reacts with Rabbit PAI-1. It binds to a conformational epitope on PAI-1 converting it into the latent form.
MOLECULAR WEIGHT:	Human PAI-1 is ~160,000
CONCENTRATION:	3.26 mg/mL
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.15M NaCl; pH 7.4
STORAGE:	Store at -70°C
IMMUNOGEN:	Human PAi-1
PURIFICATION:	Purified by immobilized Protein A.
CLASS:	IgG
SPECIES SPECIFICITY:	Human, Rabbit