



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Mouse PAI-1 (MoAB) - Capture

Product Number: PI 67

Aliquot: 1.0 mg

Lot Number 1105

Storage: -80°C

<b>DESCRIPTION:</b>	This capture monoclonal antibody produced in a PAI-1 knockout mouse.
<b>MOLECULAR WEIGHT:</b>	43 kDa non-glycosylated 50 kDa glycolysated
<b>CLONE:</b>	H34G6
<b>CONCENTRATION:</b>	1.83 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 Sodium Phosphate; 0.15 M NaCl; pH 7.4
<b>STORAGE:</b>	-80°C
<b>IMMUNOGEN:</b>	Mouse PAI-1
<b>PURIFICATION:</b>	Affinity purified by immobilized PAI-1
<b>PURITY:</b>	≥ 90% assessed by visual inspection of a Coomassie <sup>®</sup> blue-stained SDS-PAGE gel.
<b>CLASS</b>	IgG
<b>APPLICATIONS:</b>	ELISA
<b>SPECIES SPECIFICITY:</b>	Mouse