



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Human PAI-1 MoAB (non-inhibitory)

Product Number: PI 71

Aliquot: 1.0 mg

Lot Number: pi71.605A

Storage: -70°C

| | |
|-----------------------------|--|
| DESCRIPTION: | This monoclonal antibody produced in mouse and specific for human PAI-1 and exhibits non-inhibitory properties. This antibody will react with both free and complexed forms of PAI-1 |
| MOLECULAR WEIGHT: | 43 kDa non-glycosylated PAI-1 50 kDa glycosylated PAI-1 |
| CLONE: | 3A120 |
| CONCENTRATION: | 5.52 mg/mL determined spectrophotometrically. |
| STORAGE BUFFER: | 0.05 M Sodium Phosphate; 0.1 M NaCl; pH 7.4 |
| STORAGE: | -70°C. AVOID REPEATED FREEZE THAWS. |
| IMMUNOGEN: | Recombinant Human PAI-1 |
| CLASS: | IgG |
| PURIFICATION: | Purified by Protein A Sepharose. |
| APPLICATIONS: | Suitable for use in Western blotting. |
| SPECIES SPECIFICITY: | Human. |