



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Goat Anti-Rabbit PAI-1 Antiserum

Product Number: PI72

Aliquot: 1.0 mL

Lot Number: Typical Lot

Storage: -70°C

DESCRIPTION:	Anti-PAI-1 (Rabbit), Goat Polyclonal Antiserum.
MOLECULAR WEIGHT:	N/A
CONCENTRATION:	Whole antiserum
SIZE:	1.0 mL
VOLUME:	1.0 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Goat serum
APPLICATIONS:	ELISA, IHC, WB
BUFFER:	N/A

For research purposes only.