



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Protein Kinase C γ
Product Number: PK 25
Aliquot: 1.0 mL
Lot Number: 070293CT
Storage: -70°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C γ (PKC γ).
MOLECULAR WEIGHT:	78.4 kDa Calculated
CONCENTRATION:	Not Determined
ISOTYPE:	NA
STORAGE BUFFER:	Phosphate Buffered Saline; 0.05% Sodium Azide
STORAGE:	-70°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
IMMUNOGEN:	KLH conjugated synthetic peptide corresponding to amino acids 681- 697 of PKC γ . Sequence: ⁶⁸¹ V-H-P-D-A-R-S-P-I-S-P-T-P-V-P-V-M ⁶⁹⁷
PURIFICATION:	NA
APPLICATIONS:	Western Blotting: 1:100 Dilution
SPECIES SPECIFICITY:	Human, mouse and rat among other mammalian species.
REFERENCE:	JBC (1993) 268 :20110-20115