



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA 800-692-4633 • Fax 248-852-4466

www.oxfordbiomed.com

Product specifications

Anti-Protein Kinase C Delta (δ)

Product # PK 33 Typical Lot

Immunoglobulin

Type: Polyclonal Rabbit IgG

Specificity:

Prepared against synthetic peptide corresponding to a distinct C terminal region ⁷²⁶SFVNPKYEQFLE⁷³⁷ from rat protein kinase C and based on known homologies is expected to react with most mammalian species.

Applications:

Intended for use in Western blots, ELISA, Immunocytochemistry and Immunohistochemical techniques for fresh or frozen sections. Paraffin embedded or fixed tissues have not been tested.

Recommended

dilutions:

For Immunohistochemistry, immunocytochemistry and Western blotting using alkaline phosphatase or horseradish peroxidase, a dilution of 1:250 is recommended. With a more sensitive detection system, such as enhanced chemiluminescence or ELISA, the antibody may be diluted further.

Contents per vial:

100 μ L of rabbit antiserum containing 0.05% sodium azide.

Storage:

Stable for several months at 4°C. up to a year at -20°C. and indefinitely at -80° C. Avoid repeated freeze/thaw cycles.



References:

1. Wetsel, W. C., et al., *J. Cell Biol.*, **117** : 21-133 (1992)
2. *JBC*, **268** : 20110-20115.