



### Product Specification Sheet

Description: Protein Kinase C  
Epsilon Synthetic Peptide  
Product Number: PK37  
Aliquot: 0.1mg  
Lot Number: 11293JC  
Storage: -20°C

- PEPTIDE TYPE:** Distinct region of protein kinase C $\epsilon$  C-terminal  
Sequence: NH<sub>2</sub>-C-N-Q-E-E-F-K-G-F-S-Y-F-G-E-D-COOH
- CONTENTS PER VIAL:** 100  $\mu$ g of peptide in 500  $\mu$ L of ammonium bicarbonate buffer.
- STORAGE:** Stable at 4°C. For long term storage, aliquot product and freeze at -20°C. Avoid repeated freeze/thaw cycles.
- APPLICATIONS:** Will compete with the anti-peptide antibody in Western blotting for confirmation of the PKC specificity of the antibody binding.