



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product specifications

PKC Peptide Panel

Product # PK 52

This panel is an excellent way to confirm the presence of a given PKC isoform on Western blots. Each vial contains 100 µg of each of the 8 distinct C-terminal PKC peptides.

The panel contains one vial each, of the following:

PK 12 - Recombinant PKC alpha standard

PK 17 - Recombinant PKC beta I standard

PK 22 - Recombinant PKC beta II standard

PK 27 - Recombinant PKC gamma standard

PK 32 - Recombinant PKC delta standard

PK 37 - Recombinant PKC epsilon standard

PK 42 - Recombinant PKC zeta standard

PK 47 - Recombinant PKC eta standard