



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Immobilized Human Plasmin
Product Number: PL45
Aliquot: 1.0 mL resin in a 2 mL slurry.
Lot Number: pl45.1207
Storage: 4°C

DESCRIPTION:	Human Lys plasmin immobilized on agarose resin.
CONCENTRATION:	3 mg plasmin per ml of resin.
STORAGE BUFFER:	0.1M HEPES; 0.1M NaCl; 0.02% Sodium Azide; pH 7.4
STORAGE:	4°C
APPLICATIONS:	Plasmin retains catalytic activity when coupled and may be used to efficiently activate single-chain tPA and uPA to the two-chain form