



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse Antiplasmin IgG Fraction

Product Number: PL 49

Aliquot: 1.0 mg

Lot Number: pl49.506A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for mouse antiplasmin and is provided as an IgG preparation.
MOLECULAR WEIGHT:	67 kDa for mature antiplasmin**
CONCENTRATION:	2.1 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Mouse Antiplasmin
ISOTYPE:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
APPLICATIONS:	Western blotting
SPECIES SPECIFICITY:	Mouse
**REFERENCE:	<i>Hanss, Michel ML a et al., Blood Coagulation & Fibrinolysis. 14(1):107-111, January 2003.</i>