



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse Antiplasmin MoAB IgG (Inhibitory)
Product Number: PL 57
Aliquot: 1.0 mg
Lot Number: pl57.804A
Storage: -70°C

DESCRIPTION:	This mouse monoclonal antibody produced in an antiplasmin knockout mouse. The antibody is provided as an IgG prep and is inhibitory.
MOLECULAR WEIGHT:	67 kDa for mature antiplasmin**
CONCENTRATION:	1.9 mg/mL determined spectrophotometrically.
CLONE:	MAP27C9
STORAGE BUFFER:	0.5 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.8
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Mouse Antiplasmin
ISOTYPE:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
APPLICATIONS:	Western blotting
SPECIES SPECIFICITY:	Mouse
**REFERENCE:	<i>Hanss, Michel ML a et al., Blood Coagulation & Fibrinolysis.</i> 14(1):107-111, January 2003.