



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Murine Plasminogen/Plasmin Anti-Serum  
Product Number: PL 60  
Aliquot: 1.0 mg  
Lot Number: pl60.905A  
Storage: -70°C

<b>DESCRIPTION:</b>	This rabbit polyclonal antibody is generated against mouse plasminogen and will cross react with plasmin. The antibody is provided as anti-serum.
<b>CONCENTRATION:</b>	Not Applicable
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Mouse Plasminogen
<b>PURIFICATION:</b>	Not Purified
<b>APPLICATIONS:</b>	Western blotting and ELISA
<b>SPECIES SPECIFICITY:</b>	Mouse