



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse Plasminogen Affinity Purified IgG

Product Number: PL 64

Aliquot: 0.1 mg

Lot Number: pl64.105A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for mouse plasminogen and is useful in cancer and fibrinolytic studies. This is a high titer preparation.
MOLECULAR WEIGHT:	92 kDa
CONCENTRATION:	2.9 mg/mL determined spectrophotometrically.
PURITY:	> 98% pure by visual inspection on a Coomassie™ Blue stained SDS-PAGE.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 6.6
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Purified Mouse Plasminogen
CLASS:	IgG
PURIFICATION:	Affinity purified by immobilized plasmin.
APPLICATIONS:	This preparation is suitable for use in Western blotting and ELISA.
SPECIES SPECIFICITY:	Mouse