



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse Plasminogen IgG –FL Labeled

Product Number: PL 65

Aliquot: 0.5 mg

Lot Number: pl65.1104A

Storage: -70°C

DESCRIPTION:	This rabbit polyclonal antibody is specific for mouse plasminogen and is useful in cancer and fibrinolytic studies. This IgG preparation has been fluorescent labeled.
MOLECULAR WEIGHT:	92 kDa
CONCENTRATION:	3.4 mg/mL determined spectrophotometrically.
PURITY:	> 98% pure by visual inspection on a Coomassie™ Blue stained SDS-PAGE.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 6.6
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Purified Mouse Plasminogen
CLASS:	IgG
PURIFICATION:	Purified by Protein A Sepharose
APPLICATIONS:	ELISA: 1:10,000-1:100,000 dilution Western blot: 1:1,000 dilution
SPECIES SPECIFICITY:	Mouse