



OXFORD BIOMEDICAL RESEARCH
P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product Specification

Anti-Murine Plasminogen/ Plasmin (rabbit polyclonal)

Product # PL 66

Lot # asmplg-gf-bio-602

Molecular Weight: 160 KDa

Physical Appearance: colorless frozen solution

Solubility: > 10 mg/mL

Spectrophotometric Data:

Ultraviolet: Absorbance (280nm) = 2.54 0.1% = 1.4

Concentration = 8.0 mg/mL

Applications:

Western: 1:500 dilution

Considerations:

Storage Conditions: Store frozen at -70°C – stability > 1yr;

Avoid repeated freeze-thaw cycles

Buffer: 0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 7.4

Other: > 98% pure by SDS-PAGE