

Product Specification Sheet

Human Plasminogen, FITC
Labeled
Product Number: PL71
Aliquot: 1.0 mg
Lot Number: pl71.410
Storage: -70°C

DESCRIPTION: FITC Labeled Human Glu Plasminogen

MOLECULAR WEIGHT: 92,000

CONCENTRATION: 3.5 mg/mL

STORAGE BUFFER: 0.1 M HEPES, 0.1M NaCl, pH 7.4

STORAGE: -70°C Protect from light.

SOURCE: Human plasma.

LABELING: 3.0 mole FITC per mole plasminogen.

PURITY: >95% by SDS-PAGE analysis.