



Product Specification Sheet

Anti-Rat Cytochrome P450 3A1 MAb

Product Number: PM 40

Aliquot: 100 μ L

Lot Number pm40.041400sv

Storage: -80°C

DESCRIPTION:	Monoclonal Antibody Against Rat Cytochrome P450 3A1
SPECIFICITY:	Rat 3A1/3A2
APPLICATION:	Western Blotting
CONTENTS:	100 μ L per vial. Contains sufficient purified monoclonal mouse IgG for at least 10 western blots.
WORKING DILUTION:	A dilution of 1/1000 for western blots will yield a strong 3A1/3A2 band with rat liver microsomes.
STORAGE:	Store at -80°C. Avoid multiple freeze/thaw cycles.
REFERENCES:	Gelboin et al. (1985) J. biochem. Pharmacol. 34: 2225-2234 Gelboin et al. (1987) Am. Soc. Pharm. & Exp. Therap. 32: 615-624 Gelboin et al. (1993) J. biochem. Pharmacol. 46: 787-790