



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA  
USA: 800-692-4633 • Fax: 248-852-4466  
www.oxfordbiomed.com

## *Product specifications*

### **Rabbit anti-protein phosphatase 2 X/C (PP4C)**

Product # PP 17      Typical Lot

<b>Type:</b>	Polyclonal antibody
<b>Description:</b>	Purified peptide conjugated to KLH corresponding to the sequence amino acids 294 – 307.
<b>Specificity:</b>	PP X/C (34kDa)
<b>Reactivity:</b>	Human, mouse rat and bovine
<b>Presentation:</b>	100 µL whole rabbit serum provided as 0.2 µm sterile filtered solution in phosphate buffered saline with 0.08% sodium azide.
<b>Storage:</b>	Store at -20°C. Shelf life is one year. In order to avoid repeated freeze-thaw cycles, aliquot this product for individual use.
<b>Applications:</b>	Immunoblotting - Recommended dilution 1:500 - 1:1000 Immunocytochemistry - Recommended dilution 1:500 - 1:1000