



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Protein Phosphatase 2C α/β

Product Number: PP 18

Aliquot: 100 μ g

Lot Number pp18.10667A

Storage: -80°C

DESCRIPTION:	This is a rabbit polyclonal antibody is specific for serine/threonine proten phosphatase 2C α/β (PP 2C α/β) corresponding to a specific peptide immunogen. This antibody does not cross-react with other protein phosphatase isoforms.
MOLECULAR WEIGHT:	44-46 kDa
CONCENTRATION:	1 mg/mL determined by BCA method using BSA as a standard.
STORAGE BUFFER:	Phosphate Buffered Saline; 0.08% (w/v) sodium azide; Sterile filtered solution.
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	NH ₂ -Arg-Tyr-Gly-Leu-Ser-Ser-Met-Gly-Trp-Arg-Val-Glu-COOH Conjugated to KLH
CLASS:	IgG
PURIFICATION:	Ammonium Sulfate Fractionation
APPLICATIONS:	This antibody is suitable for use in Immunoprecipitation and Western blotting at 5-10 μ g/mL.
SPECIES SPECIFICITY:	Human, mouse, rat, bovine.