



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Protein Phosphatase 2 A/B γ 2

Product Number: PP 21

Aliquot: 100 μ g

Lot Number pp21.10435A

Storage: -80°C

DESCRIPTION:	This is a rabbit polyclonal antibody is specific for serine/threonine proten phosphatase 2 A/B γ 2 (PP2A/B γ 2) corresponding to amino acids 53-66 of PP2A/B. This antibody does not cross-react with other protein phosphatase isoforms.
MOLECULAR WEIGHT:	~50 kDa
CONCENTRATION:	1 mg/mL determined by BCA method using BSA as a standard.
STORAGE BUFFER:	Phosphate Buffered Saline; 0.08% (w/v) sodium azide; Sterile filtered solution.
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Amino acids 53-66 of PP2A/B
ISOTYPE:	IgG
PURIFICATION:	Ammonium Sulfate Fractionation
APPLICATIONS:	This antibody is suitable for use in Western blotting at 5-10 μ g/mL.
SPECIES SPECIFICITY:	Human, mouse, rat, bovine.