



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA
USA: 800-692-4633 • Fax: 248-852-4466
www.oxfordbiomed.com

Product specifications

Preimmune Goat Serum

Product # PR 31 C

Typical Lot

- Uses:** To be used as a negative control with immunoassays or western blots.
- Contents/vial:** 1 mL of preimmune goat serum with 0.05% sodium azide.
- Stability:** Stable for several days at 4^o C, several months at -20°C, and indefinitely at -70°C. Avoid repeated freeze/thaw cycles.