



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA: 800-692-4633 • Fax: 248-852-4466

www.oxfordbiomed.com

Product specifications

Anti-Rabbit Cytochrome P450 2E1

Product # PR 32 Typical Lot

Immunoglobulin

Type: Polyclonal goat IgG.

Specificity: Reacts as a single band with liver microsomes obtained from induced (untreated and treated with pyridine) rat or rabbit liver microsomes.

Uses: Western blotting.

Contents: Each vial contains 100 μ L of purified goat IgG. Dilution of 1/1000 for Western blots will yield a positive result with rat liver microsomes using alkaline phosphatase conjugated donkey anti-goat secondary antibodies according to the manufacturer's instructions.

Storage: Stable for several weeks at 4°C, several months at -20°C, and indefinitely at -70°C. Avoid repeated freezing and thawing.