



## Product Specification Sheet

Anti-Rabbit Cytochrome b<sub>5</sub>  
Product Number: PR91  
Aliquot: 100 µL  
Lot Number: pr91.080416  
Storage: -80° C

|                          |  |
|--------------------------|--|
| <b>IMMUNOGLOBULIN:</b>   | Polyclonal Goat IgG  |
| <b>IMMUNOGEN:</b>        | Rabbit cytochrome b <sub>5</sub>   |
| <b>CONCENTRATION:</b>    | Not determined   |
| <b>SPECIFICITY:</b>      | Specific for cytochrome b <sub>5</sub> in Rabbit, Rat, and Human; it is not cross-reactive with the 50 and 60 kDa species. |
| <b>STORAGE:</b>          | Store at -80° C. Avoid repeating freeze/thaw cycles, as it may result in loss of activity.                                 |
| <b>PURITY:</b>           | Not purified.  |
| <b>WORKING DILUTION:</b> | Suggested working dilution for Western blots is 1:5000.  |