



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Mouse tPA Antiserum

Product Number: PT 30

Aliquot: 1.0 mL

Lot Number pt30.503A

Storage: -80°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and is presented as antiserum.
<b>MOLECULAR WEIGHT:</b>	70 kDa
<b>CONCENTRATION:</b>	Not determined
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.01 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	-80°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Recombinant Mouse tPA
<b>CLASS:</b>	NA
<b>PURIFICATION:</b>	NA
<b>APPLICATIONS:</b>	This antibody exhibits utility in Western applications
<b>SPECIES SPECIFICITY:</b>	Mouse.