



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse tPA Biotin Labeled IgG

Product Number: PT 32

Aliquot: 1.0 mg

Lot Number 404

Storage: -80°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is biotin labeled.
MOLECULAR WEIGHT:	70 kDa
CONCENTRATION:	1.65 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 7.4
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Mouse tPA
CLASS:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
APPLICATIONS:	Western: 5 µg/mL
SPECIES SPECIFICITY:	Mouse.