



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Mouse tPA Affinity Purified High Titer

Product Number: PT 34

Aliquot: 0.1 mg

Lot Number pt33.904A

Storage: -80°C

| | |
|-----------------------------|---|
| DESCRIPTION: | This polyclonal antibody is produced in rabbit and has been affinity purified to produce this high titer aliquot. |
| MOLECULAR WEIGHT: | 70 kDa |
| CONCENTRATION: | 1.3 mg/mL determined spectrophotometrically. |
| STORAGE BUFFER: | 0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 7.4 |
| STORAGE: | -80°C. AVOID REPEATED FREEZE THAWS. |
| IMMUNOGEN: | Recombinant Mouse tPA |
| CLASS: | IgG |
| PURIFICATION: | Affinity purified by immobilized tPA |
| APPLICATIONS: | Western: 1:10,000 ELISA: 1:100,000 |
| SPECIES SPECIFICITY: | Mouse. |