



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Human Receptor Associate Protein MoAB

Product Number: RA 20

Aliquot: 1.0 mg

Lot Number ra20.605A

Storage: -80°C

<b>DESCRIPTION:</b>	This is a mouse monoclonal antibody is specific for the 38 kDa low density lipoprotein receptor-associated protein (RAP).
<b>MOLECULAR WEIGHT:</b>	38 kDa
<b>CLONE:</b>	MA7F1
<b>CONCENTRATION:</b>	1.08 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Tris-HCl; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	-80°C. <b>AVOID REPEATED FREEZE THAWS</b>
<b>IMMUNOGEN:</b>	Human RAP
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Purified by Protein A Sepharose.
<b>PURITY:</b>	≥ 98% assessed by visual inspection of a Coomassie <sup>®</sup> blue-stained SDS-PAGE gel.
<b>APPLICATIONS:</b>	Suitable for use in Western blot at 1 µg/mL.
<b>SPECIES SPECIFICITY:</b>	Human. Does not cross react with mouse or rabbit.