



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-85 kD LRP (Human) Biotin Labeled
Product Number: RA21
Aliquot: 0.1 mg
Lot Number: ra21.1006
Storage: -70°C

DESCRIPTION:	Biotin labeled monoclonal antibody to human 85 kD light chain of LRP.
SIZE:	0.1 mg
CONCENTRATION:	6.4 mg/mL
VOLUME:	0.016 mL
MOLECULAR WEIGHT:	160,000
BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Mouse ascites
PURIFICATION:	Purified by immobilized Protein A. IgG1 class.

For research purposes only.