

**Product Specification Sheet**

Description: Anti-Human VLDL  
Receptor Monoclonal Antibody  
Product Number: RA32  
Aliquot: 1.0 mg  
Lot Number: ra32.115  
Storage: -70°C

<b>DESCRIPTION:</b>	Anti-Human VLDL Receptor Monoclonal Antibody. Binds to an epitope on the C-terminus of VLDL Receptor (AA 191-355).
<b>MOLECULAR WEIGHT:</b>	85 kDa
<b>CONCENTRATION:</b>	4.74 mg/mL
<b>STORAGE BUFFER:</b>	0.05M Sodium Phosphate, 0.1M NaCl, 1mM EDTA, pH 6.6
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Human VLDL Receptor
<b>PURIFICATION:</b>	Immobilized protein A
<b>CLASS:</b>	IgG1
<b>APPLICATIONS:</b>	ELISA, Western Blotting (non-reducing), Blocks apoE4 binding.
<b>SPECIES SPECIFICITY:</b>	Human. Cross reactivity to other species has not been tested.