

Product Specification Sheet

Rabbit Anti-Human Prorenin, IgG Fraction
Product Number: RN26
Aliquot: 1.0 mg
Lot Number: rn26.310
Storage: -70°C

DESCRIPTION:	This polyclonal antibody was produced in rabbit and is against human prorenin. The antibody recognizes renin and prorenin.
MOLECULAR WEIGHT:	43.7 kDa
CONCENTRATION:	13.5 mg/mL
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	-70°C; minimize freeze/thaw cycles
IMMUNOGEN:	Recombinant full length wt human glycosylated protein
PURIFICATION:	IgG fraction was purified by Protein A
CLASS:	Polyclonal IgG
PURITY:	IgG Fraction
APPLICATIONS:	Western Blotting and ELISA
SPECIES SPECIFICITY:	Human