

Product Specification Sheet

Description: Rabbit anti-DMPO
Product Number: RT15
Aliquot: 100uL
Lot Number: rt15.n1331
Storage: -70°C

DESCRIPTION: Rabbit polyclonal antibody to DMPO.

IMMUNOGEN: 5, 5-dimethyl-1-pyrroline-N-oxide (DMPO)-octanoic acid adduct conjugated to ovalbumin.

PURIFICATION: Unpurified antiserum

APPLICATIONS: Recommended starting dilutions, subject to optimization:
Western Blot: 1:5,000
ELISA: 1:5,000 to 1:10,000
IHC: 1:1,000

SPECIES SPECIFICITY: None