



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Affinity purified sheep anti mouse tPA,
biotin labeled
Product Number: UP66
Aliquot: 0.1 mg
Lot Number: 607
Storage: -70°C

DESCRIPTION:	Affinity purified sheep anti mouse tPA, biotin labeled
MOLECULAR WEIGHT:	70,000
CONCENTRATION:	10.4 mg/mL
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C
IMMUNOGEN:	Recombinant Mouse tPA, full length, wild-type and glycosylated
PURIFICATION:	Affinity purified via immobilized Mouse tPA
CLASS:	Polyclonal
APPLICATIONS:	This antibody has been tested in ELISA and Western Blotting applications.
SPECIES SPECIFICITY:	Mouse