

**Product Specification Sheet**

Description: Mouse Monoclonal  
Antibody against Human uPA  
Product Number: VA54  
Aliquot: 1.0 mg  
Lot Number: VA54.typical lot  
Storage: -70°C

<b>DESCRIPTION:</b>	Mouse monoclonal antibody against the C-terminal beta (heavy) chain of human uPA. Includes the active site.
<b>MOLECULAR WEIGHT:</b>	54,000 kDa
<b>CONCENTRATION:</b>	10.0 mg/mL
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C
<b>PURIFICATION:</b>	Purified by immobilized Protein G.
<b>CLASS:</b>	IgG
<b>APPLICATIONS:</b>	ELISA, Western Blotting
<b>SPECIES SPECIFICITY:</b>	Human