

**Product Specification Sheet**

Description: Active Mouse  
Urokinase Labeled with Alexa  
Fluor 700  
Product Number: VA56  
Aliquot: 50 ug  
Lot Number: va56.912  
Storage: -70°C

<b>DESCRIPTION:</b>	Active two-chain HMW Mouse Urokinase Labeled with Alexa Fluor 700. The Urokinase was recombinantly produced in insect cells.
<b>MOLECULAR WEIGHT:</b>	44 kDa
<b>CONCENTRATION:</b>	0.49 mg/mL
<b>STORAGE BUFFER:</b>	0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA; pH 5.0
<b>STORAGE:</b>	-70°C. Protect from light.
<b>LABEL:</b>	0.74 mol Alexa Fluor 700 per 1.0 mol Urokinase
<b>PURITY:</b>	>95% by SDS-PAGE analysis